

ARMER antibody



Catalog Number: 108156

Product name

ARMER antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARMER Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, K-562 cells. Observed molecular weight by Western-blot: 22-24kd

Preparation

This antibody was obtained by immunization of ARMER recombinant protein (Accession Number: NM_015161). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:108156(ARL6IP1 antibody) at dilution of 1:300