

# ARMETL1 antibody



Catalog Number: 108157

## Product name

ARMETL1 antibody

## Specificity

Human; other species not tested.

## Antibody description

ARMETL1 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human liver cancer tissue. Positive WB detected in mouse testis tissue.

## Preparation

This antibody was obtained by immunization of ARMETL1 recombinant protein (Accession Number: NM\_001029954). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108157(ARMETL1 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:108157(ARMETL1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:108157(ARMETL1 Antibody) at dilution of 1:50 (under 40x lens)