

ARPC3 antibody



Catalog Number: 108166

Product name

Polyclonal

ARPC3 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

ARPC3 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HeLa cells, mouse brain tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 21 kDa

Dilutions

Recommended Dilution:

WB: 1:500-1:5000


Preparation

This antibody was obtained by immunization of ARPC3 recombinant protein (Accession Number: NM_001278556). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:108166(ARPC3 antibody) at dilution of 1:500

Clonality