

Arrestin C antibody



Catalog Number: 108173

Product name

Arrestin C antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Arrestin C Rabbit Polyclonal antibody. Positive IHC detected in human retinoblastoma tissue. Positive WB detected in human brain tissue, mouse eye tissue, mouse heart tissue. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of Arrestin C recombinant protein (Accession Number: BC012096). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108173(ARR3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human retinoblastoma using Catalog No:108173(ARR3 antibody) at dilution of 1:50 (under 10x lens)