

ARSG antibody

Catalog Number: 108176

Product name

ARSG antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARSG Rabbit Polyclonal antibody. Positive WB detected in human cerebellum tissue, HEK-293 cells. Observed molecular weight by Western-blot: 62 kDa

Preparation

This antibody was obtained by immunization of ARSG recombinant protein (Accession Number: XM_017024364). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:108176(ARSG antibody) at dilution of 1:300