

ARHGAP5 antibody



Catalog Number: 108178

Product name

ARHGAP5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARHGAP5 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 170-190 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001173). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108178 (ARHGAP5 antibody) at dilution of 1:800