

ARHGAP9 antibody



Catalog Number: 108179

Product name

ARHGAP9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARHGAP9 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HL-60 cells, Jurkat cells. Observed molecular weight by Western-blot: 70kd, 65kd

Preparation

This antibody was obtained by immunization of ARHGAP9 recombinant protein (Accession Number: XM_047429331). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations




Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108179(ARHGAP9 antibody) at dilution of



1:500

HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:108179(ARHGAP9 antibody) at dilution of



1:500

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108179(ARHGAP9 antibody) at dilution of 1:1000