

# ARHGEF16 antibody



Catalog Number: 108184

## Product name

ARHGEF16 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ARHGEF16 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in mouse kidney tissue. Observed molecular weight by Western-blot: 80 kDa

## Preparation

This antibody was obtained by immunization of ARHGEF16 recombinant protein (Accession Number: XM\_024446454). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108184(ARHGEF16 antibody) at dilution of



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:108184(ARHGEF16 antibody) at dilution of 1:50 (under 25x lens)