

ARL15 antibody



Catalog Number: 108197

Product name

ARL15 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARL15 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of ARL15 recombinant protein (Accession Number: NM_019087). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:108197(ARL15 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:108197(ARL15 antibody) at dilution of 1:50 (under 40x lens)