

Arylsulfatase D antibody



Catalog Number: 108208

Product name

Arylsulfatase D antibody

Specificity

Human; other species not tested.

Antibody description

Arylsulfatase D Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 65kd

Preparation

This antibody was obtained by immunization of Arylsulfatase D recombinant protein (Accession Number: BC020229). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:108208(ARSD antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:108208(ARSD antibody) at dilution of 1:50 (under 10x lens)