

ASB9 antibody



Catalog Number: 108216

Product name

Polyclonal

ASB9 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse; other species not tested.

Applications

Antibody description

ELISA, WB

ASB9 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, HepG2 cells. Observed molecular weight by Western-blot: 28k-32d

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of ASB9 recombinant protein (Accession Number: NM_001168530). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108216(ASB9 antibody) at dilution of 1:300