

# ASBT antibody

Catalog Number: 108217

## Product name

ASBT antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ASBT Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 38-40 kDa

## Preparation

This antibody was obtained by immunization of ASBT recombinant protein (Accession Number: NM\_000452). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-SLC10A2 (IP:Catalog No:108217, 4ug; Detection:Catalog No:108217 1:300) with mouse kidney tissue lysate 4000ug.