

ASCC1 antibody



Catalog Number: 108218

Product name

ASCC1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ASCC1 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, human ileum tissue. Observed molecular weight by Western-blot: 45.5 kDa

Preparation

This antibody was obtained by immunization of ASCC1 recombinant protein (Accession Number: NM_001369093). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:108218(ASCC1 Antibody) at dilution of 1:600