

# ASCC3 antibody



Catalog Number: 108220

## Product name

ASCC3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ASCC3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 251 kDa

## Preparation

This antibody was obtained by immunization of ASCC3 recombinant protein (Accession Number: NM\_022091). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108220(ASCC3 antibody) at dilution of 1:400



IP Result of anti-ASCC3 (IP:Catalog No:108220, 4ug; Detection:Catalog No:108220 1:400) with HeLa cells lysate 1080ug.