

# ARF5-Specific antibody



Catalog Number: 108248

## Product name

ARF5-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ARF5-Specific Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, COLO 320 cells, human testis tissue, mouse brain tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 21 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001662). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:108248(ARF5-Specific antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using ARF5 antibody Catalog No:108248 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).