

# ARH antibody



Catalog Number: 108257

## Product name

ARH antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ARH Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 34 kDa

## Preparation

This antibody was obtained by immunization of ARH recombinant protein (Accession Number: NM\_015627). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108257(LDLRAP1 Antibody) at dilution of 1:1000



IP Result of anti-LDLRAP1 (IP:Catalog No:108257, 4ug; Detection:Catalog No:108257 1:1000) with HEK-293 cells lysate 3000ug.