

ARHGAP20 antibody



Catalog Number: 108263

Product name

ARHGAP20 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARHGAP20 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 129-133 kDa, 88 kDa

Preparation

This antibody was obtained by immunization of ARHGAP20 recombinant protein (Accession Number: XM_005271628). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108263(ARHGAP20 antibody) at dilution of 1:500