

ASPRV1 antibody



Catalog Number: 108269

Product name

Polyclonal

ASPRV1 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

ASPRV1 Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue. Observed molecular weight by Western-blot: 28kd,37kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of ASPRV1 recombinant protein (Accession Number: BC031997). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:108269(ASPRV1 antibody) at dilution of 1:300