

# ATAD3B antibody



Catalog Number: 108278

## Product name

ATAD3B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ATAD3B Rabbit Polyclonal antibody. Positive WB detected in MAF cells. Observed molecular weight by Western-blot: 68 kDa

## Preparation

This antibody was obtained by immunization of ATAD3B recombinant protein (Accession Number: NM\_031921). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB


## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations

Mitochondria (MAFs)



WB result of ATAD3B antibody (Catalog No:108278, 1:500) with mitochondrial lysates from MAFs cells (by Maria Eugenia Soriano of UNIVERSITY OF PADUA, ITALY)