

Ataxin 2 antibody

Catalog Number: 108280

Product name

Ataxin 2 antibody

WB: 1:500-1:5000

IP: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

Ataxin 2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells, Ethacrynic acid treated HepG2 cells. Positive IHC detected in human gliomas tissue. Positive WB detected in HeLa cells, DU 145 cells, HEK-293 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 140-150 kDa

Validations



Preparation

This antibody was obtained by immunization of Ataxin 2 recombinant protein (Accession Number: NM_001310123). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108280(ATXN2 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, IHC, IF, WB, IP

Dilutions


Recommended Dilution:




Immunohistochemical of paraffin-embedded human gliomas using Catalog No:108280(ATXN2 antibody) at dilution of 1:50 (under 10x lens)

Ataxin 2 antibody


Catalog Number: 108280



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:108280(ATXN2 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of 293 cells, using ATXN2 antibody Catalog No:108280 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-Ataxin 2 (IP:Catalog No:108280, 4ug; Detection:Catalog No:108280 1:1000) with HEK-293 cells lysate 1200ug.