

ATE1 antibody



Catalog Number: 108281

Product name

ATE1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ATE1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse liver tissue. Observed molecular weight by Western-blot: 59-65kd

Preparation

This antibody was obtained by immunization of ATE1 recombinant protein (Accession Number: NM_007041). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108281(ATE1 antibody) at dilution of 1:1000