

# ATG16L1 antibody

Catalog Number: 108291

## Product name

ATG16L1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ATG16L1 Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in MCF7 cells, Jurkat cells. Observed molecular weight by Western-blot: 63-71kd

## Preparation

This antibody was obtained by immunization of ATG16L1 recombinant protein (Accession Number: BC000061). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108291(ATG16L1 antibody) at dilution of 1:1000



IP Result of anti-ATG16L1 (IP:Catalog No:108291, 4ug; Detection:Catalog No:108291 1:300) with MCF-7 cells lysate 2000ug.