

# ATG2A antibody



Catalog Number: 108293

## Product name

ATG2A antibody

## Specificity

Human; other species not tested.

## Antibody description

ATG2A Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Raji cells. Positive IHC detected in human tonsillitis tissue, human kidney tissue. Observed molecular weight by Western-blot: 215 kDa

## Preparation

This antibody was obtained by immunization of ATG2A recombinant protein (Accession Number: NM\_015104). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108293(ATG2A Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:108293(ATG2A Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108293(ATG2A Antibody) at dilution of 1:50

# ATG2A antibody



Catalog Number: 108293

---

(under 40x lens)