

ATG2A antibody

Catalog Number: 108293

Product name

ATG2A antibody

Specificity

Human; other species not tested.

Antibody description

ATG2A Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Raji cells. Positive IHC detected in human tonsillitis tissue, human kidney tissue. Observed molecular weight by Western-blot: 215 kDa

Preparation

This antibody was obtained by immunization of ATG2A recombinant protein (Accession Number: NM_015104). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108293(ATG2A Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:108293(ATG2A Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108293(ATG2A Antibody) at dilution of 1:50

ATG2A antibody



Catalog Number: 108293

(under 40x lens)