

ATG4D antibody

Catalog Number: 108298

Product name

ATG4D antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ATG4D Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 53kd

Preparation

This antibody was obtained by immunization of ATG4D recombinant protein (Accession Number: NM_032885). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108298(ATG4D antibody) at dilution of 1:800