

ATP6V1E2 antibody



Catalog Number: 108312

Product name

Polyclonal

ATP6V1E2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

ATP6V1E2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, MCF7 cells. Observed molecular weight by Western-blot: 31-33kd

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of ATP6V1E2 recombinant protein (Accession Number: XM_011533153). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108312(ATP6V1E2 antibody) at dilution of 1:400