

# ATP6V1G1 antibody



Catalog Number: 108314

## Product name

ATP6V1G1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ATP6V1G1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human kidney tissue. Positive IHC detected in human kidney. Observed molecular weight by Western-blot: 14kd

## Preparation

This antibody was obtained by immunization of ATP6V1G1 recombinant protein (Accession Number: NM\_004888). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108314(ATP6V1G1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human kidney using Catalog No:108314(ATP6V1G1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:108314(ATP6V1G1 antibody) at dilution of 1:100 (under 40x lens)