

ATP9A antibody

Catalog Number: 108318

Product name

ATP9A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ATP9A Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Positive IHC detected in human testis tissue, rat brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 105 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006045). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:108318(ATP9A antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108318(ATP9A Antibody) at dilution of 1:200 (under 10x lens).




Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108318(ATP9A Antibody) at dilution of 1:200

ATP9A antibody

Catalog Number: 108318

(under 40x lens).



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:108318(ATP9A Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)