

ATP11B antibody



Catalog Number: 108337

Product name

ATP11B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ATP11B Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 12kd

Preparation

This antibody was obtained by immunization of ATP11B recombinant protein (Accession Number: BC033880). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:108337(ATP11B antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:108337(ATP11B antibody) at dilution of 1:100 (under 40x lens)