

ATP12A antibody



Catalog Number: 108338

Product name

ATP12A antibody

Specificity

Human; other species not tested.

Antibody description

ATP12A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, A431 cells. Observed molecular weight by Western-blot: 100-110 kDa

Preparation

This antibody was obtained by immunization of ATP12A recombinant protein (Accession Number: NM_001185085). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108338(ATP12A Antibody) at dilution of 1:500