

ATP5C1 antibody



Catalog Number: 108351

Product name

ATP5C1 antibody

Immunogen

[Human ATP5C1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ATP5C1 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of ATP5C1 recombinant protein (Accession Number: XM_047425254). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:108351(ATP5C1 antibody) at dilution of 1:500