

ATP5C1 antibody



Catalog Number: 108351

Product name

ATP5C1 antibody

Clonality

Polyclonal

Immunogen

[Human ATP5C1 Recombinant protein \(GST tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

ATP5C1 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 33 kDa


Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of ATP5C1 recombinant protein (Accession Number: XM_047425254). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:108351(ATP5C1 antibody) at dilution of 1:500