

ATP5G2-Specific antibody



Catalog Number: 108354

Product name

Polyclonal

ATP5G2-Specific antibody

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, IP

Antibody description

ATP5G2-Specific Rabbit Polyclonal antibody.
Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 28 kDa

Dilutions


Recommended Dilution:

IP: 1:200-1:1000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005176).
Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

IP Result of anti-ATP5G2-Specific (IP:Catalog No:108354, 3ug; Detection:Catalog No:108354 1:300) with HEK-293 cells lysate 3200ug.