

# ATP5G2-Specific antibody

Catalog Number: 108354

## Product name

ATP5G2-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

ATP5G2-Specific Rabbit Polyclonal antibody.  
Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005176).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-ATP5G2-Specific (IP:Catalog No:108354, 3ug; Detection:Catalog No:108354 1:300) with HEK-293 cells lysate 3200ug.