

BARD1 antibody



Catalog Number: 108358

Product name

BARD1 antibody

Immunogen

Specificity

Human; other species not tested.

Antibody description

BARD1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 90-100 kDa

Preparation

This antibody was obtained by immunization of BARD1 recombinant protein (Accession Number: NM_000465). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:108358(BARD1 Antibody) at dilution of 1:300