

# BATF antibody

Catalog Number: 108362

## Product name

BATF antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BATF Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IF detected in HepG2 cells, Hela cells.

## Preparation

This antibody was obtained by immunization of BATF recombinant protein (Accession Number: NM\_006399). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of HepG2 cells, using BATF antibody Catalog No:108362 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108362(BATF Antibody) at dilution of 1:1000