

BBS2 antibody

Catalog Number: 108369

Product name

BBS2 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BBS2 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse testis tissue. Positive IF detected in hTERT-RPE1 cells. Observed molecular weight by Western-blot: 80kd

Preparation

This antibody was obtained by immunization of BBS2 recombinant protein (Accession Number: NM_031885). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108369(BBS2 antibody) at dilution of 1:500



IF result (cilia stain) of anti-BBS2 (Catalog No:108369, 1:50) with serum-starved hTERT-RPE1 cells by Dr. Moshe Kim.



IP Result of anti-BBS2 (IP:Catalog No:108369, 5ug; Detection:Catalog No:108369 1:500) with

BBS2 antibody



Catalog Number: 108369

mouse liver tissue lysate 6000ug.