

BCAS2 antibody



Catalog Number: 108380

Product name

BCAS2 antibody

Specificity

Human, Mouse, Rat, Monkey; other species not tested.

Antibody description

BCAS2 Rabbit Polyclonal antibody. Positive IF detected in Hela cells, HepG2 cells, MCF-7 cells. Positive IHC detected in human colon cancer tissue, human breast cancer tissue. Positive IP detected in mouse brain tissue. Positive WB detected in MCF7 cells, COS-7 cells, HeLa cells, hESC cells, human brain tissue, K-562 cells, mouse brain tissue, mouse liver tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 26-32kd

Preparation

This antibody was obtained by immunization of BCAS2 recombinant protein (Accession Number: NM_005872). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108380(BCAS2 antibody) at dilution of 1:1000




IP Result of anti-BCAS2 (IP:Catalog No:108380, 3ug; Detection:Catalog No:108380 1:1000) with mouse brain tissue lysate 6500ug.

BCAS2 antibody

Catalog Number: 108380



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:108380(BCAS2 antibody) at dilution of 1:100 (under 25x lens)



Immunofluorescent analysis of HeLa cells, using BCAS2 antibody Catalog No:108380 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).