

# BCAT1/ECA39 antibody



Catalog Number: 108382

## Product name

BCAT1/ECA39 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BCAT1/ECA39 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, COLO 320 cells, mouse brain tissue, mouse heart tissue. Positive IHC detected in human lymphoma tissue, human colon cancer tissue, human heart tissue, human pancreas tissue. Observed molecular weight by Western-blot: 43-45 kDa

## Preparation

This antibody was obtained by immunization of BCAT1/ECA39 recombinant protein (Accession Number: BC033864). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:108382(BCAT1 Antibody) at dilution of 1:50



(under 40x lens)

human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108382(BCAT1 antibody) at dilution of 1:300