

BCHE antibody



Catalog Number: 108385

Product name

BCHE antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

BCHE Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 85-90 kDa

Preparation

This antibody was obtained by immunization of BCHE recombinant protein (Accession Number: BC008396). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108385(BCHE Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:108385(BCHE Antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemistry of paraffin-embedded

BCHE antibody



Catalog Number: 108385

human liver tissue slide using Catalog
No:108385(BCHE Antibody) at dilution of 1:50
(under 40x lens)



Immunofluorescent analysis of HeLa cells using
Catalog No:108385(BCHE Antibody) at dilution of
1:50 and Alexa Fluor 488-conjugated AffiniPure
Goat Anti-Rabbit IgG(H+L)