

B4GALT4 antibody

Catalog Number: 108395

Product name

B4GALT4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

B4GALT4 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, human brain tissue, human kidney tissue, mouse ovary tissue. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of B4GALT4 recombinant protein (Accession Number: XM_006713799). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:108395(B4GALT4 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108395(B4GALT4 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog

B4GALT4 antibody



Catalog Number: 108395

No:108395(B4GALT4 Antibody) at dilution of 1:200 (under 40x lens).