

BACE1 antibody

Catalog Number: 108401

Product name

BACE1 antibody

Specificity

Human, Mouse, Zebrafish; other species not tested.

Antibody description

BACE1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in mouse brain tissue, HeLa cells, human brain tissue, mouse pancreas tissue, zebrafish tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 70-75 kDa

Preparation

This antibody was obtained by immunization of BACE1 recombinant protein (Accession Number: NM_012104). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200

Validations




mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108401(BACE1 antibody) at dilution of 1:500



IP Result of anti-BACE1 (IP:Catalog No:108401, 3ug; Detection:Catalog No:108401 1:500) with mouse brain tissue lysate 7500ug.

BACE1 antibody

Catalog Number: 108401



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:108401(BACE1 Antibody) at dilution of 1:50 (under 40x lens)