

# BAG3 antibody



Catalog Number: 108409

## Product name

BAG3 antibody

Recommended Dilution:

WB: 1:500-1:5000

## Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:500-1:5000


IHC: 1:20-1:200

## Antibody description

BAG3 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in HeLa cells, DU 145 cells, HEK-293 cells, K-562 cells, mouse brain tissue, mouse placenta tissue, mouse uterus tissue, rat brain tissue. Positive IHC detected in human gliomas tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 74-80 kDa

IF: 1:10-1:100

## Validations



## Preparation

This antibody was obtained by immunization of BAG3 recombinant protein (Accession Number: NM\_004281). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108409(BAG3 antibody) at dilution of 1:1000

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IF, IP, IHC


## Dilutions



IP Result of anti-BAG3 (IP:Catalog No:108409, 4ug; Detection:Catalog No:108409 1:1000) with K-562 cells lysate 11000ug.



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:108409(BAG3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:108409(BAG3 Antibody) at dilution of 1:50 (under 40x lens)

