

BAG3 antibody

Catalog Number: 108409

Product name

BAG3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BAG3 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in HeLa cells, DU 145 cells, HEK-293 cells, K-562 cells, mouse brain tissue, mouse placenta tissue, mouse uterus tissue, rat brain tissue. Positive IHC detected in human gliomas tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 74-80 kDa

Preparation

This antibody was obtained by immunization of BAG3 recombinant protein (Accession Number: NM_004281). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IP, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108409(BAG3 antibody) at dilution of 1:1000




IP Result of anti-BAG3 (IP:Catalog No:108409, 4ug; Detection:Catalog No:108409 1:1000) with K-562 cells lysate 11000ug.

BAG3 antibody


Catalog Number: 108409



Immunofluorescent analysis of HeLa cells, using BAG3 antibody Catalog No:108409 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:108409(BAG3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:108409(BAG3 Antibody) at dilution of 1:50 (under 40x lens)