

BAIAP3 antibody



Catalog Number: 108414

Product name

BAIAP3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BAIAP3 Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of BAIAP3 recombinant protein (Accession Number: NM_003933). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108414(BAIAP3 Antibody) at dilution of 1:300



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:108414(BAIAP3 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)