

# BAK antibody

Catalog Number: 108415

## Product name

BAK antibody

IHC: 1:20-1:200

IF: 1:20-1:200

## Specificity

Human, Mouse; other species not tested.

## Antibody description

BAK Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, A431 cells, RAW 264.7 cells. Positive IF detected in HeLa cells. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of BAK recombinant protein (Accession Number: NM\_001188). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IHC

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108415(BAK antibody) at dilution of 1:3000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:108415(BAK Antibody) at dilution of 1:200 (under 10x lens).




Immunohistochemistry of paraffin-embedded

# BAK antibody

Catalog Number: 108415

---

human tonsillitis tissue slide using Catalog No:108415(BAK Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of HeLa cells using Catalog No:108415(BAK Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)