

BAK antibody



Catalog Number: 108415

Product name

BAK antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

BAK Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, A431 cells, RAW 264.7 cells. Positive IF detected in HeLa cells. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of BAK recombinant protein (Accession Number: NM_001188). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108415(BAK antibody) at dilution of 1:3000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:108415(BAK Antibody) at dilution of 1:200 (under 10x lens).




Immunohistochemistry of paraffin-embedded

BAK antibody

Catalog Number: 108415

human tonsillitis tissue slide using Catalog No:108415(BAK Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of HeLa cells using Catalog No:108415(BAK Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)