

# BAMBI antibody



Catalog Number: 108416

## Product name

BAMBI antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BAMBI Rabbit Polyclonal antibody. Positive IP detected in mouse spleen tissue. Positive WB detected in rat Kidney. Observed molecular weight by Western-blot: 55kd

## Preparation

This antibody was obtained by immunization of BAMBI recombinant protein (Accession Number: NM\_012342). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IP

## Dilutions


Recommended Dilution:

WB: N/AIP: 1:200-1:1000

## Validations



WB result of Catalog No:108416 in Rat Kidney from Dr. Richard Gilbert. 55kd is Dimer of BAMBI.



IP Result of anti-BAMBI (IP:Catalog No:108416, 3ug; Detection:Catalog No:108416 1:300) with mouse spleen tissue lysate 4000ug.