

# BAP1 antibody



Catalog Number: 108419

---

## Product name

BAP1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

BAP1 Rabbit Polyclonal antibody. Positive WB detected in A2780 cells, human placenta tissue, mouse ovary tissue. Observed molecular weight by Western-blot: 91 kDa

## Preparation

This antibody was obtained by immunization of BAP1 recombinant protein (Accession Number: NM\_004656). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:108419(BAP1 antibody) at dilution of 1:300