

BCL2L10 antibody



Catalog Number: 108423

Product name

BCL2L10 antibody

Specificity

Human; other species not tested.

Antibody description

BCL2L10 Rabbit Polyclonal antibody. Positive WB detected in A2780 cells. Positive IHC detected in human lung tissue, human kidney tissue. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of BCL2L10 recombinant protein (Accession Number: NM_020396). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:117096(BCL2L10 antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human lung using Catalog No:117096(BCL2L10 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung using Catalog No:117096(BCL2L10 antibody) at dilution of 1:50 (under 40x lens)